

STERILE SURGICAL BLADES & SCALPELS




SURGICAL KNIFE HANDLES
Plastic - Autoclavable - Reusable



STERILE DISPOSABLE SURGICAL BLADES





Kai surgical blades have a complete selection of the most widely used blade size - both in **carbon steel** and **stainless steel**. Sterility is guaranteed by use of gamma radiation. 100 blades individually foil packed, hermetically sealed, are packed in a dispenser box.



STERILE DISPOSABLE SCALPELS

- Stainless steel blade
- With transparent cap
- Sterilized by E.O.G.
- 20 scalpels in a dispenser box



SURGICAL KNIFE HANDLES

Autoclavable-Reusable-Plastic

- Polypropylene plastic body
- With easy blade remover
- Each packed in polybag
- 10 handles in a box

Handle type: Style NO.3 P-3


Handle type: Style NO.4 P-4



Product specifications may be changed without notice

DISPOSABLE SAFETY SCALPELS

STERILE • STAINLESS STEEL • SINGLE USE



Features

- Safe from accidental needle sticks
- Blade extends and retracts with just one hand
- Safely returns exposed blade to original covered position
- Securely locks blade into both extended and retracted position
- Glides easily and comfortably into locked position
- Metric ruler assists in determining length of incision
- Sterile and ready to use
- Disposable

Package

- 10 scalpels in a dispenser box
- 100 boxes per shipping carton

To extend blade



To retract blade




There are two sizes of handles depending on the blade shapes.

Small : SS-10, SS-11, SS15

Large : SS-12, SS-20, SS-21, SS-22, SS-23, SS-24, SS-25

SS-10




SS-11



SS-14




SS-12




SS-15




SS-20




SS-21



SS-24




SS-22



SS-25



SS-23



Product specifications may be changed without notice.



Certified Management System
ISO 9001
ISO 14001
ISO 13485



kai industries co.,ltd.

1110,Oyana, Seki City, Gifu Pref. 501-3992 Japan
Phone:+81 (0) 575 28 6600 Fax:+81 (0) 575 28 6611
Email:info-ind@kai-group.com <http://www.kai-ind.co.jp>



kai Europe GmbH

Löhdorfer Str. 51, D-42699 Solingen, Germany
Phone:+49 (0) 212 232380 Fax:+49 (0) 212 2323899

PRINTED IN JAPAN KIM-011 0211 / 0510 K 2K1